

derma+flex[®]

gel tissue adhesive

DermaFlex is a patented medical adhesive formulated as a controlled flow gel.

Dermaflex is sterile. It is not runny. It goes and stays where you want it to.

While multiple layers of adhesive can be applied it is not recommended as DermaFlex's higher viscosity allows one layer to be applied. The thicker the layer of adhesive the longer the curing time.


Dermaflex contains no solvents, will not sting and its unique polymerization method (curing) reduces the risk of an exothermic reaction (heat generation) which is considered to be negligible.

As the name suggests DermaFlex provides a secure, invisible and flexible seal that acts as a protectant over most closed incisions, minor wounds and abrasions.



DERMA+ FLEX DAUBER TIP

Allows a flow controlled application to the surface, that provides a flexible waterproof bandage



DERMA+ FLEX FINE APPLICATOR

Allows fine applications where looks (cosmetic issues) are considered

DermaFlex is ideally suited for areas of the skin exposed to flexing, moisture and difficult to seal area's with traditional bandages such as between fingers and toes, fingertips, chins, eyebrows, elbows and knees.

Dermaflex comes with two applicators. A nozzle for fine applications where looks (cosmetic issues) are considered and a unique patented dauber top for larger area's and control of the adhesive protectant layer.

Chemence Medical's "Sure+Close" wound system augments the applications for this unique Product.

For more information or on how to order please contact:
www.MedicalGlue.Com



Waterproof & Sting Free

